



MIQE checklist for authors, reviewers and editors

ITEM TO CHECK	IMPORTANCE	CHECKLIST
EXPERIMENTAL DESIGN		
Definition of experimental and control groups	E	
Number within each group	E	
Assay carried out by core lab or investigator's lab?	D	
Acknowledgement of authors' contributions	D	
SAMPLE		
Description	E	
Volume/mass of sample processed	D	
Microdissection or macrodissection	E	
Processing procedure	E	
If frozen - how and how quickly?	E	
If fixed - with what, how quickly?	E	
Sample storage conditions and duration (especially for FFPE samples)	E	
NUCLEIC ACID EXTRACTION		
Procedure and/or instrumentation	E	
Name of kit and details of any modifications	E	
Source of additional reagents used	D	
Details of DNase or RNase treatment	E	
Contamination assessment (DNA or RNA)	E	
Nucleic acid quantification	E	
Instrument and method	E	
Purity (A260/A280)	D	
Yield	D	
RNA integrity method/instrument	E	
RIN/RQI or Cq of 3' and 5' transcripts	E	
Electrophoresis traces	D	
Inhibition testing (Cq dilutions, spike or other)	E	
REVERSE TRANSCRIPTION		
Complete reaction conditions	E	
Amount of RNA and reaction volume	E	
Priming oligonucleotide (if using GSP) and concentration	E	
Reverse transcriptase and concentration	E	
Temperature and time	E	
Manufacturer of reagents and catalogue numbers	D	
Cqs with and without RT	D*	
Storage conditions of cDNA	D	
qPCR TARGET INFORMATION		
If multiplex, efficiency and LOD of each assay.	E	
Sequence accession number	E	
Location of amplicon	D	
Amplicon length	E	





qPCR VALIDATION		
Evidence of optimisation (from gradients)	D	
Specificity (gel, sequence, melt, or digest)	E	
For SYBR Green I, Cq of the NTC	E	
Standard curves with slope and y-intercept	E	
PCR efficiency calculated from slope	E	
Confidence interval for PCR efficiency or standard error	D	
r2 of standard curve	E	
Linear dynamic range	E	
Cq variation at lower limit	E	
Confidence intervals throughout range	D	
Evidence for limit of detection	E	
If multiplex, efficiency and LOD of each assay.	E	
DATA ANALYSIS		
qPCR analysis program (source, version)	E	
Cq method determination	E	
Outlier identification and disposition	E	
Results of NTCs	E	
Justification of number and choice of reference genes	E	
Description of normalisation method	E	
Number and concordance of biological replicates	D	
Number and stage (RT or qPCR) of technical replicates	E	
Repeatability (intra-assay variation)	E	
Reproducibility (inter-assay variation, %CV)	D	
Power analysis	D	
Statistical methods for result significance	E	
Software (source, version)	E	
Cq or raw data submission using RDML	D	

Table 1. MIQE checklist for authors, reviewers and editors. All essential information (E) must be submitted with the manuscript. Desirable information (D) should be submitted if available. If using primers obtained from RTPrimerDB, information on qPCR target, oligonucleotides, protocols and validation is available from that source.

*: Assessing the absence of DNA using a no RT assay is essential when first extracting RNA. Once the sample has been validated as RDNA-free, inclusion of a no-RT control is desirable, but no longer essential.

** : Disclosure of the probe sequence is highly desirable and strongly encouraged. However, since not all commercial pre-designed assay vendors provide this information, it cannot be an essential requirement. Use of such assays is advised against.

