

TripleA Pipette Tip Compatibility Chart



TripleA Pipette Tips are certified to be consistent in form and function. This chart is based on internal research and is subject to change. Therefore it is recommended to test the compatibility of TripleA tips on your pipettes.

				5010, 5110, 5210, 5510, 5610	5020, 5030	5040, 5140, 5240, 5340, 5540, 5640	5050, 5150, 5155, 5250, 5550, 5650	5060, 5360	5070, 5170, 5175, 5270, 5570, 5670
Brand	Model	Volume Range	Channels	10µl	20µl, 100µl	200µl	300µl	1 ml	1250µl
Biohit	mLine	0.1-3µl, 0.5-10µl	1	•					
Biohit	mLine	2-20µl, 10-100µl, 20-200µl	1		•	•	•		
Biohit	mLine	100-1000µl	1					•	•
Biohit	mLine	0.5-10µl	8/12	•					
Biohit	mLine	5-100µl	8/12		•	•	•		
Biohit	mLine	30-300µl	8/12				•		
Biohit	Proline	0.1-2.5µl, 0.5-10µl	1	•					
Biohit	Proline	2-20µl, 5-50µl, 10-100µl, 20-200µl	1		•	•	•		
Biohit	Proline	100-1000µl	1					•	•
Biohit	Proline	0.5-10µl	8/12	•					
Biohit	Proline	5-50µl	8/12		•	•	•		
Biohit	Proline	50-300µl	8/12				•		
Brand	Transferpette	0.1-1µl, 0.5-10µl	1	•					
Brand	Transferpette	2-20µl micro	1	•					



				5010, 5110, 5210, 5510, 5610	5020, 5030	5040, 5140, 5240, 5340, 5540, 5640	5050, 5150, 5155, 5250, 5550, 5650	5060, 5360	5070, 5170, 5175, 5270, 5570, 5670
Brand	Model	Volume Range	Channels	10µl	20µl, 100µl	200µl	300µl	1 ml	1250µl
Brand	Transferpette	2-20µl, 5-50µl, 10-100µl, 20-200µl	1		•	•	•		
Brand	Transferpette	100-1000µl	1					•	•
Brand	Transferpette	0.5-10µl	8/12	•					
Brand	Transferpette	2-20µl	8/12	•					
Brand	Transferpette	2.5-25µl, 5-50µl, 10-100µl, 20-200µl	8/12		•	•	•		
Brand	Transferpette	30-300µl	8/12				•		
Brand	Transferpette S	0.1-1µl, 0.1-2.5µl, 0.5-10µl	1	•					
Brand	Transferpette S	2-20µl, 10-100µl, 20-200µl	1		•	•	•		
Brand	Transferpette S	100-1000µl	1					•	•
Brand	Transferpette S	0.5-10µl	8/12	•					
Brand	Transferpette S	5-50µl, 10-100µl, 20-200µl	8/12		•	•	•		
Brand	Transferpette S	30-300µl	8/12				•		
Eppendorf	Research Plus	0.1-2.5µl, 0.5-10µl	1	•					
Eppendorf	Research Plus	2-20µl micro	1	•					
Eppendorf	Research Plus	2-20µl, 10-100µl, 20-200µl	1		•	•	•		
Eppendorf	Research Plus	30-300µl	1				•		
Eppendorf	Research Plus	100-1000µl	1					•	•
Eppendorf	Research Plus	0.5-10µl	8/12	•					
Eppendorf	Research Plus	10-100µl	8/12		•	•	•		
Eppendorf	Research Plus	30-300µl	8/12				•		
Eppendorf	Reference	0.1-2.5µl, 0.5-10µl	1	•					
Eppendorf	Reference	2-20µl micro	1	•					
Eppendorf	Reference	2-20µl, 10-100µl, 50-200µl	1		•	•	•		
Eppendorf	Reference	100-1000µl	1					•	•



				5010, 5110, 5210, 5510, 5610	5020, 5030	5040, 5140, 5240, 5340, 5540, 5640	5050, 5150, 5155, 5250, 5550, 5650	5060, 5360	5070, 5170, 5175, 5270, 5570, 5670
Brand	Model	Volume Range	Channels	10µl	20µl, 100µl	200µl	300µl	1 ml	1250µl
Thermo	Finnpette 4500	0.2-2µl, micro, 0.5-5µl, micro	1	•					
Thermo	Finnpette 4500	0.5-5µl, 2-20µl, 5-50µl, 10-100µl, 20-200µl	1		•	•	•		
Thermo	Finnpette 4500	100-1000µl	1					•	•
Thermo	Finnpette 4500	0.5-10µl	8/12	•					
Thermo	Finnpette 4500	5-50µl	8/12		•	•	•		
Thermo	Finnpette 4500	50-300µl	8/12				•		
Gilson	Pipetman P2 Classic/G/L	0.2-2ul, 1-10ul	1	•					
Gilson	Pipetman P20 Classic/G/L	2-20ul, 20-100ul, 50-200ul	1		•	•	•		
Gilson	Pipetman P1000 Classic/G/L	200-1000ul	1					•	•
Gilson	Pipetman P10L	0.5-10µl, 2-20µl, 20-200µl	8/12		•	•	•		
Gilson	Pipetman P300L	20-300µl	8/12				•		
Nichiryo	Nichipet Premium/EX/EX Plus	0.1-2µl	1	•					
Nichiryo	Nichipet Premium/EX/EX Plus	2-20µl, 10-100µl, 20-200µl	1		•	•	•		
Nichiryo	Nichipet Premium/EX/EX Plus	100-1000µl	1					•	•
Nichiryo	Nichipet 7000	0.5-10µl	8/12	•					
Nichiryo	Nichipet 7000	5-50µl, 40-200µl	8/12		•	•	•		
Nichiryo	Nichipet 7000	50-300µl	8/12				•		
Rainin	Pipet-Lite SL-2XLS	0.1-2µl, 0.5-10µl	1	•					
Rainin	Pipet-Lite SL-20XLS	2-20µl, 10-100µl, 20-200µl	1		•	•	•		
Rainin	Pipet-Lite SL-300XLS	20-300µl	1				•		
Rainin	Pipet-Lite SL-1000XLS	100-1000µl	1					•	•
Rainin	Classic PR-2	0.1-2µl, 0.5-10µl	1	•					
Rainin	Classic PR-20	2-20µl, 10-100µl, 20-200µl	1		•	•	•		



				5010, 5110, 5210, 5510, 5610	5020, 5030	5040, 5140, 5240, 5340, 5540, 5640	5050, 5150, 5155, 5250, 5550, 5650	5060, 5360	5070, 5170, 5175, 5270, 5570, 5670
Brand	Model	Volume Range	Channels	10µl	20µl, 100µl	200µl	300µl	1 ml	1250µl
Rainin	Classic PR-1000	100-1000µl	1					•	•
Socorex	Calibra 822	0.2-2µl, 1-10µl micro	1	•					
Socorex	Calibra 822	1-10µl, 2-20µl, 10-100µl, 20-200µl	1		•	•	•		
Socorex	Calibra 822	100-1000µl	1					•	•
Socorex	Calibra 852	1-10µl	8	•					
Socorex	Calibra 852	10-100µl, 20-200µl	8/12		•	•	•		
Socorex	Acura 825	0.1-2µl, 0.5-10µl	1	•					
Socorex	Acura 825	1-10µl, 2-20µl, 5-50µl, 10-100µl, 20-200µl	1		•	•	•		
Socorex	Acura 825	100-1000µl	1					•	•
Socorex	Acura 855	0.5-10µl	8/12	•					
Socorex	Acura 855	5-50µl, 20-200µl	8/12		•	•	•		
Socorex	Acura 855	40-350µl	8/12				•		

